



## RHINOCELL® EPISTAXIS PACKING

### Description

RHINOCELL® Epistaxis Packs are constructed of a unique sponge material using an exclusive formula designed for rapid and effective treatment of epistaxis. RHINOCELL® Epistaxis Packs provide anterior epistaxis control (EP-551525) and posterior epistaxis control (EP-1001525). The small cell structure of RHINOCELL® provides a smooth surface, which maximizes comfort while minimizing tissue ingrowth. The high tensile strength of the sponge material prevents shedding or tearing. The packs are supplied compressed for ease of insertion and will expand as fluid is introduced, giving time for accurate positioning. All RHINOCELL® Epistaxis Packs are supplied with an integral lanyard to secure dressing to the outside nose.

### Indications for use

For use as nasal epistaxis packing.

### Supplied

Sterile

10 packs per dispenser carton, individually pouched



EP-551525



EP-1001525

### FEATURES

- Biocompatible lint and fiber-free sponge material
- Superior liquid absorption and wicking
- Supplied compressed for ease of insertion
- Extremely soft and pliable cell structure
- Attached lanyard to secure packing

### Product Codes & Dimensions

REF	RhinoCell® Epistaxis Pack	
EP-551525	5.5 x 1.5 x 2.5 cm	anterior
EP-1001525	10.0 x 1.5 x 2.5 cm	posterior